

3mm Cannulae

For endometrial biopsy

EFFECTIVE · DISPOSABLE · AFFORDABLE



3mm Cannulae

- Disposable and sterile packed
- Permanently integrated winged base
- Two opposing apertures
- Approximately 23 cm (9 in) in length
- Depth markings printed in cm on cannulae
- Compatible with Ipas MVA Plus® Aspirator with a 6mm adaptor
- Compatible with Ipas MVA Single Valve Aspirator with no adaptor required
- Medical grade, latex-free copolyester plastic

Clinical Indications

- Endometrial Biopsy

Clinical Applications

- Evaluation of infertility
- Evaluation of amenorrhea
- Evaluation of suspected endometrial cancer
- Evaluation of abnormal uterine bleeding
- Diagnosis of endometrial infections

In less than 2% of cases one or more of the following can occur during or after the procedure: uterine or cervical injury, pelvic infection, vagal reaction, incomplete evacuation, or acute hematometra

**ALSO AVAILABLE FROM
WOMANCARE GLOBAL**



Single Valve Aspirator

For uterine evacuation and endometrial biopsy

- Easy disassembly and reassembly
- 60 mL volume
- 610-660 mm (24-26 in) Hg vacuum
- Latex-free plastic

Product Code: SVSS10 (sterile pack of 10)

