

Flexible Karman Cannulae

The original flexible cannulae

EFFECTIVE · DISPOSABLE · AFFORDABLE



Flexible Karman Cannulae

- Flexible
- Disposable and sterile packed
- Available in eight diameters
- Color-coded by size (adaptor kit required)
- Depth markings in cm printed on cannulae
- Two opposing apertures for cannulae with diameters of 4-8 mm
- Large, single scoop aperture for cannulae with diameters of 9, 10, and 12 mm
- Approximately 24 cm (9.5 in) in length
- Medical grade, latex-free polypropylene plastic
- FDA-cleared and compliant with ISO 13485 and CE-marked

Product Code/Diameter Chart

Size	4 mm	5 mm	6 mm	7 mm	8 mm	9 mm	10 mm	12 mm
Code	FC4	FC5	FC6	FC7	FC8	FC9	FC10	FC12



ALSO AVAILABLE FROM WOMANCARE GLOBAL



Cannulae Adaptor Sets

For use in conjunction with applicable cannulae and MVA Aspirators
 Product Codes: ADPSET (sizes 6-10 MM)
 ADP6MM-ADP10MM (set of 5)



Ipas MVA Plus® Aspirator

For uterine evacuation
 Product Code: PLUS



Denniston Dilators

For cervical dilation associated with uterine evacuation
 Product Code: DD5

THE IPAS ASPIRATOR PRODUCT LINE

A COST EFFECTIVE SOLUTION FOR UTERINE EVACUATION¹

- Backed by 40 years of global market use and Ipas' long-standing quality reputation
- FDA-cleared, CE-marked, and ISO 13485 compliant
- Small, inexpensive, portable, quiet, and safe^{2,3,4}
- Ideal for performing procedures in the outpatient setting
- Ideal for settings with no electricity or with space constraints
- Limited investment required
- Easy to use

DEMONSTRATED CLINICAL BENEFITS IN UTERINE EVACUATION

	Ipas Aspirators	EVA	Sharp Curettage
Safe	X	X	
Lower Incidence of Side Effects/Complications	X	X	
Cost Effective	X		X*
Portable	X		X
Quiet	X		X
Single Use	X		
Reusable	X	X	X
Outpatient/No Anesthesia	X	X	
Quick Recovery Time	X	X	
No Electricity Required	X		X
Requires Limited Space	X		X
Requires Limited Investment	X		X

*When used without anesthesia

PRODUCT PORTFOLIO

ASPIRATORS

- Ipas Plus[®] Aspirator
- Single Valve Aspirator
- Double Valve Aspirator

CANNULAE

- Ipas Easy Grip[®] Cannulae
- Flexible Karman Cannulae
- Rigid Uterine Cannulae

DILATORS AND ACCESSORIES

- Adapter sets, dilators and aspirator accessories also available

KITS

- Uterine Evacuation Kit

TRAINING TOOLS

- Pelvic Model

References: 1. Veruy DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial. *S Afr Med J* 83:13-15. 2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. *Int J Gynecol Obstet*, 1994, 45: 261-267. 3. De Jonge ETM, Pattinson, RC, Makin ID et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. *S Afr Med J*, 84:481-483. 4. Fawcus S, McIntyre J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. *S Afr Med J*, 87:438-442. 5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition, 2012. Annex 5, Recommendation 2, page 113. 6. Hale RW, Kobara TV, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

WomanCare Global International is registered as a charity in England and Wales (No.1138896) and a company limited by guarantee (No.07279512)
WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920)



WomanCare Global
Quality • Innovation • Choice

Phone (call free on): 011 44 20 7921 0555
Email: customerservice@womancareglobal.org
www.womancareglobal.org