

# Rigid Uterine Cannulae

The control of a sharp curette combined with the safety of manual uterine aspiration

STRAIGHT & CURVED · DISPOSABLE · AFFORDABLE



## Rigid Uterine Cannulae

- Available in both straight and curved
- For use with either Ipas aspirators or tubing for electric suction machines
- Disposable and sterile packed
- Available in 7 sizes color-coded by diameter: 6, 7, 8, 9, 10, 11 and 12 mm
- Medical grade, latex-free styrene plastic
- FDA-cleared, CE-marked, and ISO 13485 compliant

## Clinical Indications

- Treatment of incomplete abortion for uterine sizes up to 12 weeks since the woman's last menstrual period



Product Code/Size Chart

Diameter	Product Code	
	Straight	Curved
6 mm	RS6	RC6
7 mm	RS7	RC7
8 mm	RS8	RC8
9 mm	RS9	RC9
10 mm	RS10	RC10
11 mm	RS11	RC11
12 mm	RS12	RC12

## ALSO AVAILABLE FROM WOMANCARE GLOBAL



### Ipas MVA Plus<sup>®</sup> Aspirator

For uterine evacuation  
Product Code: PLUS



### Collection Tubing Set

For use with electric suction machines  
Product Code: TS1

# THE IPAS ASPIRATOR PRODUCT LINE

A COST EFFECTIVE SOLUTION FOR UTERINE EVACUATION<sup>1</sup>

- Backed by 40 years of global market use and Ipas' long-standing quality reputation
- FDA-cleared, CE-marked, and ISO 13485 compliant
- Small, inexpensive, portable, quiet, and safe<sup>2,3,4</sup>
- Ideal for performing procedures in the outpatient setting
- Ideal for settings with no electricity or with space constraints
- Limited investment required
- Easy to use

## DEMONSTRATED CLINICAL BENEFITS IN UTERINE EVACUATION

	Ipas Aspirators	EVA	Sharp Curettage
Safe	X	X	
Lower Incidence of Side Effects/Complications	X	X	
Cost Effective	X		X*
Portable	X		X
Quiet	X		X
Single Use	X		
Reusable	X	X	X
Outpatient/No Anesthesia	X	X	
Quick Recovery Time	X	X	
No Electricity Required	X		X
Requires Limited Space	X		X
Requires Limited Investment	X		X

\*When used without anesthesia

## PRODUCT PORTFOLIO

### ASPIRATORS

- Ipas Plus<sup>®</sup> Aspirator
- Single Valve Aspirator
- Double Valve Aspirator

### CANNULAE

- Ipas Easy Grip<sup>®</sup> Cannulae
- Flexible Karman Cannulae
- Rigid Uterine Cannulae

### DILATORS AND ACCESSORIES

- Adapter sets, dilators and aspirator accessories also available

### KITS

- Uterine Evacuation Kit

### TRAINING TOOLS

- Pelvic Model

References: 1. Veruy DAA, Crowther CA. (1993). Suction v. conventional curettage in complete abortion: a randomized clinical trial. *S Afr Med J* 83:13-15. 2. Blumenthal, Paul D and R.E. Remsburg. 1994. A time and cost analysis of the management of incomplete abortion with manual vacuum aspiration. *Int J Gynecol Obstet*, 1994, 45: 261-267. 3. De Jonge ETM, Pattinson, RC, Makin ID et al. (1994). Is ward evacuation for uncomplicated incomplete abortion under systematic analgesia safe and effective: a randomized clinical trial. *S Afr Med J*, 84:481-483. 4. Fawcus S, McIntyre J, Jewkes RK et al. (1997). Management of incomplete abortions at South African public hospitals. *S Afr Med J*, 87:438-442. 5. World Health Organization, Department of Reproductive Health and Research. Safe Abortion: Technical and policy guidance for health systems. Second edition. 2012. Annex 5, Recommendation 2, page 113. 6. Hale RW, Kobara TV, Sharma SD, Tsuei JJ, Gramlich EP, Nakayama RT. Office termination

WomanCare Global International is registered as a charity in England and Wales (No.1138896) and a company limited by guarantee (No.07279512)  
WomanCare Global Trading CIC is a UK Registered Community Interest Company (No.07231920)



**WomanCare Global**  
Quality • Innovation • Choice

Phone (call free on): 011 44 20 7921 0555  
Email: [customerservice@womancareglobal.org](mailto:customerservice@womancareglobal.org)  
[www.womancareglobal.org](http://www.womancareglobal.org)